

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

REF 370, **LOT** 40360 Level 1 - 3 📅 2021-02-28

REF		Units Einheiten Unitades Unités	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
HUMAN Hormones & Tumor Markers								
			Level 1 - 40361		Level 2 - 40362		Level 3 - 40363	
52010	AFP	ng/ml	19.2	12.9 - 25.5	98.8	63.1 - 134	221	148 - 296
52020	CEA	ng/ml	1.90	1.28 - 2.53	12.8	8.57 - 17.0	21.1	13.7 - 28.5
55050	Cortisol	ng/ml	37.3	23.9 - 50.7	213	137 - 291	355	227 - 483
55060	DHEA-S	µg/ml	0.52	0.29 - 0.75	1.32	0.72 - 1.91	4.68	2.57 - 6.78
55030	Estradiol	pg/ml	183	103 - 264	545	305 - 785	944	528 - 1359
55040	Estriol	ng/ml	1.81	1.00 - 2.63	7.10	3.90 - 10.3	14.6	8.00 - 21.1
52035	fPSA	ng/ml	0.42	0.28 - 0.56	1.80	1.21 - 2.39	*	-
53020	FSH	IU/l	8.66	5.80 - 11.52	20.4	13.7 - 27.1	47.1	31.6 - 63.7
54015	ft3	pg/ml	2.70	1.81 - 3.59	5.17	3.46 - 6.87	7.33	4.91 - 9.75
54025	ft4	ng/dl	0.81	0.54 - 1.07	2.29	1.54 - 3.05	2.91	1.95 - 3.86
53040	hCG	IU/l	5.21	2.89 - 7.53	16.9	11.3 - 22.4	> 193	136 - > 250
53010	LH	IU/l	1.56	1.05 - 2.08	26.13	17.5 - 34.8	84.8	56.8 - 113
53030	Progesterone	ng/ml	1.03	0.56 - 1.49	10.9	5.97 - 15.7	> 36.4	20.0 - >40.0
55020	Prolactin	ng/ml	8.01	5.37 - 10.7	22.7	15.2 - 30.2	47.1	31.5 - 62.6
52030	PSA	ng/ml	1.02	0.56 - 1.48	4.84	3.15 - 6.54	24.5	15.9 - 33.0
54010	T3	ng/ml	0.87	0.58 - 1.15	2.24	1.49 - 2.99	3.77	2.52 - 5.01
54020	T4	µg/dl	5.43	3.64 - 7.22	11.2	7.48 - 14.9	13.1	8.78 - 17.4
55010	Testosterone	ng/ml	2.01	1.17 - 2.86	4.49	2.60 - 6.37	9.41	5.46 - 13.4
54030	TSH	mIU/l	0.31	0.17 - 0.44	4.81	3.13 - 6.49	*	-

* The marked level is not applicable for the indicated parameter as the target value is out of the measuring range of the assay.

